10-minute IHC for Frozen Tissue: Step by Step Procedure

1) Fixation

- After frozen tissue is cut onto slides, immediately fix for **1 minute** in Novodix ihc Fixative or reagent grade acetone
(Alternative fixative solutions can be used based on user preference)
- Place slides to soak in wash buffer until ready to proceed

2) Blocker

- 1st wash (wash buffer) – Use spray bottle to wash tissue on the slides
- Wipe off excess wash buffer around the tissue with paper towel or Kim Wipe
- Apply blocker (covering entire tissue) and incubate for **1 minute**
DO NOT WASH OFF BLOCKER
- Hold slide vertical on a hard surface and tap off blocker onto a paper towel
- Wipe off excess blocker around the tissue with paper towel or Kim Wipe

3) Antibody

- Shake/mix antibody thoroughly before use
- Apply antibody (covering entire tissue) and incubate for **3 minutes**
- 2nd wash (wash buffer) – Use spray bottle to wash tissue on the slides
- Wipe off excess wash buffer around the tissue with paper towel or Kim Wipe

4) DAB

- Prepare DAB working solution (1 drop Chromogen to 1 ml Diluent) using provided supplies: **Before opening DAB chromogen shake or tap down fluid to the bottom of vial to ensure that it is not in the cap.**
- Use a plastic transfer pipette to dispense 1 ml of Diluent into an Eppendorf tube
- Use a plastic transfer pipette to dispense 1 drop (equal to 30 µl) of Chromogen
- Mix/invert well before use
- Apply DAB working solution (covering entire tissue) and incubate for **3 minutes**
- 3rd wash (DI or distilled water) -- Use spray bottle to wash tissue on the slides
- Wipe off excess water around the tissue with paper towel or Kim Wipe

5) Hematoxylin and Coverslip

- Apply Hematoxylin (by drop or soak) and incubate for a time range of **15 seconds up to 1 minute** (time depends on strength of user's Hematoxylin dye)
- 4th wash (DI or distilled water) – Use spray bottle to wash off Hematoxylin
- Wipe off excess water around the tissue with paper towel or Kim Wipe
- Apply mounting medium (aqueous or permanent) and Coverslip of user preference

Novodix Bulk Reagent Formulations:

1. ihc Fixative, (375ml of methyl alcohol, 100ml of 37% formaldehyde and 25ml of glacial acetic acid).

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